






THUASNE

Rebel Reliever®



fr	Genouillère rigide de décharge	4
en	Offloading rigid knee brace.....	5
de	Starre Kniestütze	6
nl	Rigide ontlastende kniebrace	8
it	Ginocchiera di scarico rigida	9
es	Rodillera rígida de descarga.....	10
pt	Joelheira rígida de alívio do peso	11
da	Stift, aflastende knæbind	13
fi	Kuormitusta keventävä jäykkä polvituki.....	14
sv	Styvt knästöd för avlastning.....	15
el	Άκαμπτη επιγονατίδα αποφόρτισης	16
cs	Pevná odlehčující kolenní ortéza.....	17
pl	Sztywna, odciążająca orteza stawu kolanowego	19
lv	Stingra atslogojoša ceļa ortoze.....	20
lt	Standus apkrovos mažinimo antkelis.....	21
et	Koormust vähendav elastne põlve ortoos.....	22
sl	Toga razbremenilna opornica za koleno	23
sk	Pevná odľahčovacia kolenná ortéza	25
hu	Merev tehermentesítő térdrögzítő	26
bg	Твърда наколенка за разтоварване.....	27
ro	Genunchieră rigidă pentru descărcare.....	28
ru	Жесткий наколенник для снятия напряжения	29
hr	Fiksna ortoza, potpora za koljena.....	31
zh	減压用硬质护膝	32
ar	داعم ركبة صلب لتخفيف الثقل.....	34


	8	9
fr	Extension	Flexion
en	Extension	Flexion
de	Extension	Flexion
nl	Extensie	Flexie
it	Estensione	Flessione
es	Extensión	Flexión
pt	Extensão	Flexão
da	Ekstension	Fleksion
fi	Ekstensio	Fleksio
sv	Extension	Flexion
el	Εκταση	Κάμψη
cs	Extenze	Flexe
pl	Wyprost	Zgięcie
lv	Atliekšana	Saliekšana
it	Istiesimas	Sulenkimas
et	Sirutus	Painutus
sl	Iztegovanje	Upogibanje
sk	Extenzia	Flexia
hu	Nyújtás	Hajlítás
bg	Разгъване	Сгъване
ro	Extensi3n	Flexi3n
ru	Разгибание	Сгибание
hr	Ekstenzija	Fleksija
zh	伸展	屈曲
ar	المَدّ	النتي

Fitting videos (English version)

fr	Décharge	Contrôle du mouvement
en	Off-loading	Motion control
de	Entlastung	Bewegungsführung
nl	Losmaken	Bewegingscontrole
it	Scarico	Controllo del movimento
es	Descarga	Control del movimiento
pt	Alívio do peso	Controlo do movimento
da	Aflastning	Kontrol af bevægelsen
fi	Vapautus	Liikkeen hallinta
sv	Avlastning	Rörelsekontroll
el	Αποφόρτιση	Έλεγχος της κίνησης
cs	Snížení zátěže	Kontrola pohybu
pl	Odciażenie	Kontrola ruchu
lv	Atslogošana	Kustības kontrole
lt	Apkrovis sumažinimas	Judesio kontrolė
et	Koormuse vähendamine	Liikuvuse kontroll
sl	Razbremenitev	Nadzor gibanja
sk	Odlahčenie	Vymedzenie rozsahu pohybu
hu	Tehermentesítés	Mozgásszabályozás
bg	Отбремявяване	Контрол на движението
ro	Descărcare	Controlul mişcării
ru	Разгрузка	Контроль движения
hr	Rasterećenje	Kontrola pokreta
zh	减压	运动控制
ar	تخفيف نقل	التحكم في الحركة

	(a)		(b)		(c)	
	cm	in	cm	in	cm	in
XS	31 - 39	12 ½ - 15 ½	7,5 - 9	3 - 3 ½	28 - 32	11 - 12 ½
SM	39 - 47	15 ½ - 18 ½	9 - 10	3 ½ - 4	31 - 35	12 ¼ - 13 ¼
MD	47 - 53	18 ½ - 21	10 - 11,5	4 - 4 ½	34 - 38	13 ¼ - 15
LG	53 - 60	21 - 23 ½	11,5 - 12,5	4 ½ - 5	37 - 40	14 ¼ - 15 ¾
XL	60 - 64	23 ½ - 25	12,5 - 14	5 - 5 ½	38 - 43	15 - 17
XXL	64 - 71	25 - 28	14 - 15,5	5 ½ - 6	43 - 48	17 - 19

fr	(a) Circonférence de la cuisse (15 cm au-dessus du genou)	(b) Largeur du genou	(c) Circonférence du mollet (15 cm au-dessous du genou)
en	(a) Thigh circumference (15 cm/6" above the knee)	(b) Knee width	(c) Calf circumference (15 cm/6" below the knee)
de	(a) Oberschenkelumfang (15 cm über dem Knie)	(b) Kniebreite	(c) Wadenumfang (15 cm unter dem Knie)
nl	(a) Omvang van de dij (15 cm boven de knie)	(b) Breedte van de knie	(c) Omvang van de kuit (15 cm onder de knie)
it	(a) Circonferenza della coscia (15 cm al di sopra del ginocchio)	(b) Larghezza del ginocchio	(c) Circonferenza del polpaccio (15 cm al di sotto del ginocchio)
es	(a) Circunferencia del muslo (15 cm por encima de la rodilla)	(b) Ancho de la rodilla	(c) Circunferencia de la pantorrilla (15 cm por debajo de la rodilla)
pt	(a) Circunferência da coxa (15 cm acima do joelho)	(b) Largura do joelho	(c) Circunferência da barriga da perna (15 cm abaixo do joelho)
da	(a) Lårets omkreds (15 cm over knæet)	(b) Knæets bredde	(c) Læggenes omkreds (15 cm under knæet)
fi	(a) Reiden ympärys (15 cm polven yläpuolelta)	(b) Polven leveys	(c) Pohkeen ympärys (15 cm polven alapuolelta)
sv	(a) Omkrets runt låret (15 cm ovanför knäet)	(b) Knäets bredd	(c) Omkrets runt vad (15 cm under knäet)
el	(a) Περίμετρος του μηρού (15 εκ. πάνω από το γόνατο)	(b) Περίμετρος του γόνατος	(c) Περίμετρος της κνήμης (15 εκ. κάτω από το γόνατο)
cs	(a) Obvod stehna (15 cm nad kolenem)	(b) Obvod kolena	(c) Obvod lýtky (15 cm pod kolenem)
pl	(a) Obwód uda (15 cm powyżej kolana)	(b) Szerokość kolana	(c) Obwód łydki (15 cm poniżej kolana)
lv	(a) Augšstilba apkārtmērs (15 cm virs ceļgala)	(b) Ceļgala platums	(c) Apakšstilba apkārtmērs (15 cm zem ceļgala)
lt	(a) Šlaunies apimtis (15 cm virš kelio)	(b) Kelio plotis	(c) Blauzdos apimtis (15 cm žemiau kelio)
et	(a) Reie ümbermõõt (15 cm põlvest kõrgemal)	(b) Põlve ümbermõõt	(c) Sääremarja ümbermõõt (15 cm põlvest allpool)
sl	(a) Obseg stegna (15 cm nad kolenom)	(b) Širina kolena	(c) Obseg goleni (15 cm pod kolenom)
sk	(a) Obvod stehna (15 cm nad kolenom)	(b) Šírka kolena	(c) Obvod lýtky (15 cm pod kolenom)
hu	(a) Comb körfogata (15 cm a térd felett)	(b) Térd szélessége	(c) Vádli körfogata (15 cm a térd alatt)
bg	(a) Обиколка на бедрото (15 cm над коляното)	(b) Ширина на коляното	(c) Обиколка на прасеца (15 cm под коляното)
ro	(a) Circumferința coapsei (la 15 cm deasupra genunchiului)	(b) Lățimea genunchiului	(c) Circumferința gambei (la 15 cm sub genunchi)
ru	(a) Диаметр окружности бедра (на 15 см выше колена)	(b) Ширина колена	(c) Обхват голени (на 15 см ниже колена)
hr	(a) Opseg bedra (15 cm iznad koljena)	(b) Širina koljena	(c) Opseg lista (15 cm ispod koljena)
zh	(a) 大腿周长 (膝盖上方15厘米处)	(b) 膝关节宽度	(c) 小腿周长 (膝盖下方15厘米处)
ar	(a) محيط الفخذ (15 سم فوق الركبة)	(b) عرض الركبة	(c) محيط ربة الساق (15 سم فوق الركبة)

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Description/Destination

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OFFLOADING RIGID KNEE BRACE

Description/Destination

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table.

Product available in left or right version.

Available in 6 standard sizes and 9 adapted sizes.

An "adapted" size is when the calf and thigh circumferences do not correspond in the standard size chart. In which case, it is possible to use a product with a different thigh and calf size.

The device is composed of:

- a frame ensuring the brace's rigidity,
- the TM5+ hinge ① mimicking the knee's natural movement, protected by a cover,
- a Loadshifter mechanism ② (tool-less) to extend the brace uprights,
- 2 condylar pads,
- 4 posterior half-straps (Ⓐ, Ⓞ, Ⓜ, Ⓠ) (including a synergistic suspension strap Ⓞ), and 2 anterior half-strap (Ⓢ, Ⓡ) ①
- 4 quick buckles to ease opening and closing of the straps (Ⓐ, Ⓞ, Ⓜ, Ⓠ, Ⓢ, Ⓡ) ②

The flexion/extension adjustment system is designed to be used without tools. Additional condylar pads, extension & flexion limitations and extra-self-fastening attachments are included in a plastic bag in the brace's box. Depending on country of sale, additional accessories/spare parts could be available.

Composition

Textile components: elastane - polyamide - ethylene vinyl acetate - silicone - poly-methyl methacrylate - polyethylene.
Rigid components: aluminium - stainless steel - polyoxymethylene - vinyl - polypropylene - polyamide - epoxy resin - brass.

Properties/Mode of action

The Rebel Reliever® knee brace is primarily designed to apply corrective forces and alignment support to your leg by means of a three-point pressure offloading system.

These corrective forces help reduce varus or valgus and help distribute excessive load away from the damaged knee compartment. Through correction of alignment, the brace will help reduce load, aiding in the reduction of compression and inflammation.

Reducing the load on the affected side would lead to the reduction of pain and cartilage degeneration.

The corrective properties of the Rebel Reliever® brace are designed to help with varus or valgus correction (+/-9°).

In neutral position, without angle correction, the Rebel Reliever® brace can also be used for a specific immobilisation in flexion or extension for the conservative treatment of knee ligament injuries or post-operative immobilisation.

- Extension can be adjusted to 0°, 5°, 10°, 15°, 20°, 30° and 40°.
- Flexion can be adjusted to 0°, 30°, 45°, 60°, 75° and 90°.

The brace maintains proper fit on the leg through:

- the Synergistic Suspension Strap Ⓞ,
- the shape of the tibial bolster (rounded lateral part, flat medial part) that prevents rotation,
- the Compression/Suspension pack providing improved hold on the thigh (optional).

Indications

Symptomatic unicompartmental femorotibial osteoarthritis (moderate to severe).
Knee off-loading for post-traumatic, post-operative or degenerative conditions.
Joint instability/laxity.

Alternative to osteotomy or leg misalignment surgery.
Conservative treatment of knee ligament injuries and/or ruptures (cruciate and/or lateral ligaments).

Contraindications

- Do not use the product if the diagnosis has not been confirmed.
- Do not use in the event of known allergy to any of the components.
- Do not apply the product in direct contact with broken skin or an open wound without an adequate dressing.
- Do not use in the event of severe varicose veins (preventing the regular use of an offloading brace).
- Do not use in the event of significant genu recurvatum.
- Do not use in the event of significant genu varum or genu valgum.
- Do not use in the event of major venous thromboembolic history without thromboprophylaxis.

Precautions

- Verify the product's integrity before every use.
- Do not use the product if it is damaged.
- Choose the appropriate size to fit the patient, referring to the size chart.
- It is recommended that a healthcare professional supervises the first application.
- Strictly comply with your healthcare professional's prescription and recommendations for use.

This product is intended for the treatment of a given condition. Its duration of use is limited to this treatment only.

For hygiene and performance reasons, do not re-use the product for another patient.

It is recommended to adequately tighten the device to achieve support/immobilisation without restricting blood circulation.

In the event of history of venous or lymphatic disorders, consult a healthcare professional.

In the event of discomfort, significant hindrance, pain, variation in limb volume, abnormal sensations or change in colour of the extremities, remove the device and consult a healthcare professional.

In the event of any modification in the product's performance, remove it and consult a healthcare professional.

Before any sports activity, check the compatibility of the use of this medical device with your healthcare professional.

Do not wear the product in a medical imaging machine.

Do not use the device in case of application of certain products on the skin (creams, ointments, oils, gels, patches...).

Undesirable side-effect

This device can cause skin reactions (redness, itching, burns, blisters, etc.) or wounds of various degrees of severity.

Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Instructions for use/Application

It is recommended to wear the device directly on the skin, unless contraindicated.

HEALTHCARE PROFESSIONAL INSTRUCTION GUIDE:

Please read these instructions with the patient when fitting the brace for the first time and make sure that he/she understands how to position the knee brace.

How to fit the knee brace:

- **Step 1: Evaluation of initial pain**
Assess the level of pain felt by the patient to measure the function and efficacy of the knee brace after it is fitted.
- **Step 2: Positioning the brace on the leg**
- Check that the 4 buckles on the back straps are open ①.
Confirm the Loadshifter Mechanism above each knee joint are in the factory-set neutral position (same height on both sides of the thigh shell).
- Have the patient sit on the edge of a chair with their knee bent to approximately 30-60° of flexion.
- Place the knee brace on the bare leg ②.
Ensure that the brace size is appropriate by checking that there is no space between the brace and the leg.
If there is any gapping, or if the brace is too tight, if possible, conform the frame of the brace or retake measurements of the leg and select the appropriate size (see size table).
- **Step 3: Fitting the knee brace**
- Close the suspension strap ④ with the buckle so it is above the calf muscle, just below the flexion fold on the back of knee ④.
- Then adjust the tightening with the self-fastening on the strap.
- Next fasten the anterior tibial strap ⑧ in a similar way, but not too much.
The anterior-posterior position of the hinges can be modulated by adjusting straps ② and ⑧.
The hinges and uprights should be aligned with the midline of the sides of the leg, or just posterior to midline (2nd third of the leg in the sagittal plane).
After you pre-set the anterior strap, instruct the patient to leave it closed when putting on or taking off the brace.
- Close the lower strap ⑥, then the two rear thigh straps ② and ⑧ and finally close the front thigh strap ⑧ with their respective buckles and self-fastenings.
Strap chafes can be angled to ensure the straps fit flat against the back side of the leg.
If one or more of the knee brace straps is too long, each strap can be cut to the desired length.
To do this, remove the self-fastening hook tab from the end of the strap, cut the strap to the desired length, and reattach the self-fastening hook tab to the end of the strap ④.
Take care to not cut any strap too short.
The comfort pads fixed inside the straps may need to be removed for this operation and repositioned after the strap is cut to avoid any interference when the strap is tightened.
- **Step 4: Check the fit, suspension and pain level prior to adjusting the corrective force**
Ask the patient to stand and walk, taking normal steps and looking straight ahead.
Confirm the brace fits and suspends properly on the leg.
Ask the patient if the pain level with the brace on, still in the neutral position, is the same or less than without the brace.
Once the brace has been adjusted and the straps properly fitted, leave the self-fastenings in position and use the buckles to place and remove the product.
During the initial fitting, there is no need to change the positioning of the self-fastening used to attach the straps to the rigid frame of the brace. They are intended for strap changes when necessary.
This operation must be done by a healthcare professional.
Refer to the instruction of use of the corresponding replacement kit.
- **Step 5: Adjusting the corrective force**
The dual Loadshifter Mechanism ② enables you to change the angle of

the femoral shell of the brace to increase the corrective 3-point force and offload the compromised (degenerated, compressed) medial or lateral compartment of the knee. ④

Both uprights ② are set during fabrication in a neutral position at mid-height.

The best way to adjust the correction is to increase the length of the upright on the affected side. This produces more leverage force.

If the patient is short, you can also increase correction by shortening the length of the upright on the unaffected side.

For severe OA and/or severe varus or valgus deformities, you may need to adjust both uprights.

To adjust the Loadshifter ③, lift the plastic Snap Lock lever located above the hinge, on the side requiring adjustment, lengthen or shorten the upright and close the lever. ②

 **Tip: if opening of the Snap Lock seems difficult, do a quick back and forth movement with both hands on the upright to unlock the mechanism. ②**

Medial Compartment Femoro-Tibial Osteoarthritis of the Knee: this brace is effective for the treatment of unicompartiment internal/medial femoral-tibial osteoarthritis of the knee.

A thicker condylar pad should be positioned on the external/lateral hinge (on the side opposite the affected compartment).

To increase corrective force, the internal/medial upright (affected side) can be lengthened compared with the external/lateral upright, or the external/lateral upright can be shortened compared to the internal/medial upright.

Lateral Compartment Femoro-Tibial Osteoarthritis of the Knee: this brace is also effective for the treatment of unicompartiment external/lateral compartment femoral-tibial osteoarthritis of the knee.

The thicker condylar pad should be positioned on the internal/medial hinge (on the side opposite the affected compartment). To increase corrective force, the internal/medial upright (affected side) can be lengthened compared with the external/lateral upright, or the external/lateral upright can be shortened compared to the internal/medial upright.

•Step 6: Evaluation of pain after adjusting the corrective force

After adjusting the offloading of the brace:

- Ask the patient to walk again and assess his/her pain level.
- Repeat the operation, increasing or reducing the correction until the proper amount of offloading is achieved for the patient.
- If the patient feels any discomfort, reduce the angle of the thigh shell which will reduce the corrective force.
- Once the adjustment of the knee brace is finalized, it is possible to trim any excess padding at the lower end of the pads, just above both hinges.

For conservative treatment, start with a small amount of correction and let the patient wear the knee brace for one or two weeks.

If it is necessary to increase the correction, see the patient again.

A bag containing an additional set of condylar pads is provided in the box with the brace.

If necessary, use the set of thicker pads to increase the compression on the sides of the knee.

Lines and numbers are printed on the uprights of the knee brace. There is no correlation between these graduations and specific degrees of correction. You can record in the patient's chart the initial setting of the Loadshifter.

Flexion/extension control:

The brace's default extension setting is 0°.

To change this setting, please follow the following instructions, to be repeated with the same limitation on both hinges.

The extension and flexion limitations are situated on plastic pieces contained in the box ④ ⑤.

Extension can be limited to 0°, 5°, 10°, 15°, 20°, 30° and 40°.

Flexion can be limited to 0°, 30°, 45°, 60°, 75° and 90°.

•Adjusting the extension limitation:

1. Choose the desired extension limitation on the corresponding plastic piece ④.
 2. Open the hinge cap by using the latch located at the back: push the latch inwards and then lift the cap.
 3. Flex the hinge slightly and remove the existing extension limitation. Take care to note the direction this limitation is facing.
 4. Insert the new extension limitation. Take care to check that the shape of this part follows the contours of the cap. Straighten the hinge to full extension to confirm the extension limitation is correctly positioned.
 5. Close the cap.
- Do a few flexions/extensions to make sure that the limitation is properly locked at the desired angle.

•Adjusting the flexion limitation:

1. Choose the desired flexion limitation on the corresponding plastic piece ⑤.
 2. Open the hinge cap by using the latch located at the back: push the latch inwards and then lift the cap.
 3. Straighten the hinge to full extension and insert the desired flexion limitation (the flexion limitations should be placed behind the hinge). Ensure the hole in the flexion limitation is aligned with the hole in the hinge (metal part).
 4. Close the cap.
- Do a few flexions/extensions to make sure that the limitation is properly locked at the desired angle.

•Definitive locking of the hinge cap (to prevent from being opened by the patient):

1. Remove from its support, the plastic piece which locks the cap ⑥.
2. When the cap is closed, position the locking piece so that it fits into the cap notch. Then push the locking piece until it comes into abutment.

Warning:

The adjustment of the flexion/extension must be decided and performed by the healthcare professional, not by the patient.

Both hinges MUST be adjusted to the same angle. There can be damage to the hinges, and compromises (including injury) to the patient, if the settings are not the same on both hinges.

PATIENT INSTRUCTIONS GUIDE:

The healthcare professional who performed the initial fitting of your brace should have made all the necessary modifications to the corrective force applied by the knee brace, selected the appropriate hinge pads, and adjusted the length of the straps. The healthcare professional should have also explained how to put on your brace.

Fitting the orthosis:

- Step 1: Sit on the edge of a chair and slightly bend your leg (30° to 60°).
- Step 2: Open the 4 buckles on the posterior straps (④, ⑤, ⑥, ⑦) and place the brace on the leg. ②

The condylar pads attached to the inside of the hinges should press against the sides of your knee.

The middle of the pads should line up between the middle and the upper third of your knee cap (bone at the front of the knee joint) and should also be centered on the side of your leg ②.

NOTE: It is always better to put your brace on a little too high than too low.

- Step 3: First close the Synergy Suspension Strap ④ above the calf using the buckle ①.

This strap should be tightened in the flexion fold behind the knee.

It is the most important strap for preventing the brace from slipping down your leg.

By tightening the synergistic suspension strap ④ above the calf muscle in the flexion fold behind the knee you will ensure that the hinges are at the correct level on the leg.

- Step 4: Next fasten and tighten the lower strap ③, then the two straps behind the thigh ② and ①.

If you feel the need to tighten the straps during an activity, make sure you tighten them following the steps indicated in the fitting instructions. It should not normally be necessary to readjust the front straps ③ and ②. Do not change the position of the self-gripping clamps used to attach the straps to the rigid frame of the brace. They are intended for strap changes when necessary. This operation must be performed by a healthcare professional.

It may take several weeks to feel comfortable with the brace on your leg. It is recommended to initially wear the brace for only a few hours a day.

Adjusting the correction:

During the initial fitting, the healthcare professional should make all necessary adjustments to obtain a base level of corrective force, according to the medical indication.

If you feel the need for additional correction in the future, or if you are feeling any discomfort, or are experiencing any other issue with the fit, function or suspension of the brace, stop wearing your brace until you can see your physician or the healthcare professional who fit the brace. Thuasne cannot be held responsible for undesirable effects or injuries resulting from any unsupervised or inadapted adjustments.

Care/Maintenance

Close the self-fastening tabs before washing. Hand wash. Do not use detergents, fabric softeners or aggressive products (products containing chlorine). Squeeze out excess water. Dry away from any direct heat source (radiator, sun, etc.). If the device is exposed to seawater or chlorinated water, make sure to rinse it in clear water and dry it.

Hinges ②:

The hinges on the brace are pre-lubricated in the factory. If sand, dirt or water gets inside the hinges, they may require lubricating again. If you notice the hinges not gliding smoothly, a few drops of a synthetic lubricant can be applied. Wipe off any excess lubricant before wearing the brace to prevent stains on clothing.

Straps:

After prolonged use, if the fibers on your strap do not adhere as well to the self-fastening tab, cut the strap shorter so the self-fastening tab adheres to a section of the strap that has fresher fibers. If this is not possible, you should contact the medical provider who fit your brace.

Pads:

The brace is lined with padding that provides a comfortable interface between the leg and the shells. The straps also have pads. Do not remove these pads from the brace or straps. Wipe the pads after each use to remove any moisture and let the pads air dry. You can also clean the pads with a mild anti-bacterial soap and rinse them off with fresh water. Do not wash the pads in a machine and do not dry in a tumble dryer.

Storage

Store at room temperature, preferably in the original packaging.

Disposal

Dispose of in accordance with local regulations.

COMMERCIAL WARRANTY AGREEMENT AND WARRANTY LIMITATIONS

Thuasne offers a free, limited commercial warranty to the user, in the territory where the device was purchased, against defects in manufacturing and workmanship for a period of:

- six months for the textile components;
- one year for the rigid components.

The limited warranty is effective from the date of purchase of the product by the end-user.

The limited commercial warranty does not apply to any defects in manufacturing and workmanship in case of:

- misuse of the product or any damage occurred by a usage outside the normal and intended use of the product as mentioned in the instructions for use;
- damages occurred as part of attempts to modify the product.

Any deterioration or improper trimming of the product during its modification or adjustment process by the healthcare professional when fitting device is expressly excluded from this warranty.

Any claim for this commercial warranty must be sent by the user to the entity where the product was purchased, which will forward this claim to the corresponding Thuasne entity.

Any warranty claim will first be reviewed by Thuasne to determine if the conditions of the limited warranty are fulfilled and do not fall into one of the cases of exclusion of the commercial warranty.

To benefit from the warranty, the buyer must mandatorily provide an original and dated proof of purchase of the product.

If the conditions of the limited warranty are fulfilled and the claim is made by the user or its legal representative (parents, guardian...) within the warranty delays indicated above, the buyer will get a new substitution product.

It is expressly agreed that this commercial warranty is in addition to the legal warranties binding the entity which sold the product to the user, in

accordance with the applicable local legislation in the country of purchase of the product.

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• Étape 4 : Fermer et serrer ensuite la sangle postérieure inférieure ©, puis les deux sangles postérieures de la cuisse ® et ⑥.
THUASNE UK Ltd